

Actolind® w Solution


Actolind® w Solution is a cleansing, irrigation, moisturising, and care solution used together with materials such as bandages, dressings, pads and wound fillers for acute, chronic and infected wounds as well as 1st and 2nd degree burns.

We are next to you in each step of wound management


We offer optimal solutions for wound healing **thanks to modern wound care dressings** designed with the clinical experience we have accumulated in wound care.




Actolind® Hydrocolloid
Hydrocolloid adhesive,
absorbent, waterproof
wound dressing



**Actolind® Hydrocolloid
Foam**
Highly retentive hydrocolloid backed
foam dressing



Actolind® Foam
Adhesive / Non Adhesive
Highly retentive foam dressing



Actolind® Wound Pad
Adhesive / Transparent
Highly retentive wound dressing

Actolind® w Solution

Article Number	Unit	Packaging Type	Packaging Unit	PZN
06.6703.05	10x5 ml	Monodose Strip	100 x 5	16889754
06.6703.10	5x10 ml	Monodose Strip	50 x 5	16684727
06.6703.50	50 ml	Bottle with spray cap	100	10985385
06.6703.100	100 ml	Bottle with spray cap	81	10985391
06.6703.250	250 ml	Bottle with spray cap	49	10985416
06.6703.350	350 ml	Bottle	30	16684733
06.6703.500	500 ml	Bottle	30	10985422
06.6703.1	1000 ml	Bottle	10	10985439



ACTO GmbH & ACTOLIND GmbH
Schulstrasse 8 • 31613 Wietzen-Holte / Germany
Tel: +49 (0) 5022 9801-0 • Fax: +49 (0) 5022 9801-37
Braunschweig Office
Büchnerstrasse 11 • 38118 Braunschweig / Germany
Bayern Branch
Balthasar - Neumann Str. 4 / 97204 Höchberg / Germany
www.actogmbh.com • info@actogmbh.com



(28.09.2021) REV_02_EN


Actolind® w Solution

Antimicrobial irrigation solution for wound cleansing,
moisturisation and care

Polyhexanide + **Poloxamer**

It accelerates the healing of **wounds and burns**

✓ **Cleansing** + ✓ **Moisturising** + ✓ **Antimicrobial Effect**



German Company




Actolind® w Solution

Actolind® w Solution is a irrigation, cleansing, moisturising and care product for acute, chronic and infected wounds as well as 1st and 2nd degree burns.

The ideal partner to **support** an easy and pain-free wound care



Effective

Cleanses the wound and the peri-wound

- Removes contaminants, scabs, necrotic tissues and biofilm layers from wounds effectively by decreasing the surface tension
- Ensures cleaning, moisturising and infection protection when applied before dressing change



Patient-friendly

Convenient and painless dressing changes

- Softens encrusted dressing materials that are dried and adhered to the wound
- Minimizes adhesion of dressings to the wound surface
- Allows atraumatic change of dressing materials



Protective

Provides antimicrobial cleansing

- Prevents possible re-infections
- Prevents the formation of biofilms
- Highly effective against multi-resistant organisms (**MRSA, VRE**) and mycotic pathogens (*C. albicans*)



Healing

Provides moisture and care

- Keeps the wound surface clean and moist
- Enhances the effectiveness of the dressing
- Prevents the adhesion of microorganisms
- Accelerates the healing of wounds and burns



Supportive

Efficient wound care thanks to the synergistic effect of poloxamer and polyhexanide (PHMB)



Dermatologically tested

Neither irritating nor painful

- Mucous membrane compatible
- Does not cause sensitization



Gentle

High biocompatibility

- Neither inhibits granulation nor epithelialization
- High tissue compatibility
- Suitable for long-term treatment
- Well tolerated by allergic and diabetic patients



Indications

- Infected or uninfected acute wounds
- Infected or uninfected chronic wounds
- Infected or uninfected thermal or chemical wounds (e.g. 1st and 2nd degree burns, chemical burns, etc.)
- Wounds caused by electrical current
- Radiation wounds
- Entry regions of urological catheters
- Maintenance of peristomal skin, fistulas and abscesses

Application Recommendations

- **Actolind® w Solution** should be used for wound surface cleansing and irrigation before any treatment.
- **Actolind® w Solution** can be heated to body temperature prior to use.
- Bandages, dressings, pads, etc. such as **Actolind® Sterile Gauze Compress**, soaked in **Actolind® w Solution** can be used to clean wounds.
- For optimal wound cleansing, it is recommended to wet the infected wound consisting of deep and intensively necrotic tissues, for at least 10-15 min. with a fully wetted compress

Compatible
with all modern
wound dressings