

WUNDVERBÄNDE UND EIGENSCHAFTEN

Produktmerkmale	Actiind® Hydrocolloid / Thin	Actiind® Hydrocolloid Foam	Actiind® Hydrocolloid Plus	Actiind® Plus Hydrocolloid	Actiind® Foam Non Adhesive	Actiind® Foam Adhesive	Actiind® Plus Foam	Actiind® Plus Foam Adhesive	Actiind® Wound Foam / Thin	Actiind® Wound Pad /Thin	Actiind® Wound Pad	Actiind® Wound Pad - Adhesive	Actiind® Wound Transparent Adhesive	Actiind® Sterile Gauze Compress	Actiind® Sterile Bandage
CE 93/42/EWG Klassifizierung	IIB	IIB	IIB	IIB	IIB	IIB	IIB	IIB	IIB	IIB	IIB	IS	IS	IIA	IS
Einmalgebrauch, steril	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Kleberand			✓			✓						✓	✓		
Leichte Anpassung an die Wunde	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Saugfähigkeit	✓ ❖	✓ ❖	✓ ❖	✓ ❖	✓ ❖	✓ ❖	✓ ❖	✓ ❖	✓ ❖	✓ ❖				✓ ❖ ❖	
Zum Zuschneiden geeignet	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Durch Kontakt mit Exsudat Gelbfärbung	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Wasserdichte, atmungsaktive Schicht	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Einzelteil-Ablösung	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Selbstklebend	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Verschiedene Größen	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Dünne Option	✓									✓					
Barrierefunktion											✓	✓	✓	✓	✓

❖ Hydrokolloid: Niedrig-Mittel
❖ Foam: Multi-Medium
❖ Gazekompress: Niedrig